

Batch Number: **B 555**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 3 December 2015
Expiration Date: 1 April 2016
Volume: 0.578 ml +/- 6 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	1.0
DAPI staining:		
Cryptosporidium % +ve	96	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 20 October 2015
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

