

Batch Number: **B 575**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 18 August 2016
Expiration Date: 16 December 2016
Volume: 0.596 ml +/- 19 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

	Mean	St.Dev.
Cryptosporidium count:	100	1.2

DAPI staining:

Cryptosporidium % +ve	100 %
------------------------------	--------------

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 30 June 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

