

Batch Number: **B 599**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 25 May 2017  
**Expiration Date:** 22 September 2017  
**Volume:** 0.598 ml +/- 16 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

**Counts**

(Method Ref: ESCG-001)

	Mean	St.Dev.
<b>Cryptosporidium count:</b>	<b>99</b>	<b>1.3</b>

**DAPI staining:**

<b>Cryptosporidium % +ve</b>	<b>100 %</b>
------------------------------	--------------

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 14 April 2017  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Officer:**



.....  
Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

