

Batch Number: **B 628**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 1 March 2018  
**Expiration Date:** 29 June 2018  
**Volume:** 0.607 ml +/- 18 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

| <b>Counts</b>                 |            |            |
|-------------------------------|------------|------------|
| <i>(Method Ref: ESCG-001)</i> | Mean       | St.Dev.    |
| <b>Cryptosporidium count:</b> | <b>100</b> | <b>1.5</b> |
| <b>DAPI staining:</b>         |            |            |
| <b>Cryptosporidium % +ve</b>  | <b>100</b> | <b>%</b>   |

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 15 February 2018  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Officer:**

.....  
Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer  
Accreditation Number: 14993  
Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

