

Batch Number: B 642

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 16 August 2018
Expiration Date: 14 December 2018
Volume: 0.579 ml +/- 1 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	1.8
DAPI staining:		
Cryptosporidium % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 17 July 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer



.....
Madhuli Patil
 28 August 2018



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:
 BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9859
www.btfbio.com

