

Batch Number: B 646

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 13 September 2018
Expiration Date: 11 January 2019
Volume: 0.618 ml +/- 18 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

	Mean	St.Dev.
Cryptosporidium count:	99	2.1

DAPI staining:

Cryptosporidium % +ve	100	%
------------------------------	------------	----------

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 25 August 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Specialist



Linh Le

25 September 2018



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9859
www.btfbio.com

