

**Batch Number: B 661**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 31 January 2019  
**Expiration Date:** 31 May 2019  
**Volume:** 0.553 ml +/- 61 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

<b>Counts</b>		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
<b>Cryptosporidium count:</b>	<b>99</b>	<b>2.1</b>
<b>DAPI staining:</b>		
<b>Cryptosporidium % +ve</b>	<b>100</b>	<b>%</b>

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 12 January 2019  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Manager**



Peter Verwoerd  
 14 February 2019



Accredited for compliance with ISO 17034  
 Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquiries:**  
 BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8221 9859  
[www.btfbio.com](http://www.btfbio.com)

