

Batch Number: B 663

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 28 February 2019
Expiration Date: 28 June 2019
Volume: 0.610 ml +/- 14 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

| Counts | | |
|-------------------------------|------------|------------|
| <i>(Method Ref: ESCG-001)</i> | Mean | St.Dev. |
| Cryptosporidium count: | 99 | 1.8 |
| DAPI staining: | | |
| Cryptosporidium % +ve | 100 | % |

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 09 January 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Specialist

Linh Le
15 March 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquiries:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9859
www.btfbio.com

