

**Batch Number: B 668**

## General

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 26 April 2019  
**Expiration Date:** 24 August 2019  
**Volume:** 0.571 ml +/- 13 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

<b>Counts</b>		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
<b>Cryptosporidium count:</b>	<b>99</b>	<b>1.7</b>
<b>DAPI staining:</b>		
<b>Cryptosporidium % +ve</b>	<b>100</b>	<b>%</b>

## Stock specifics

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 27 March 2019  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

## Quality Assurance Specialist



Linh Le  
10 May 2019



Accredited for compliance with ISO 17034  
 Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



## Water Industry Enquiries:

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8221 9859  
[www.btfbio.com](http://www.btfbio.com)

