



EasyStain™

- EasyStain™ is a kit for fluorescently staining *Cryptosporidium* and *Giardia* (C&G)
- Kit contains 4mL working strength EasyStain™ solution (80 tests), fixing buffer, positive control and slide mounting medium
- Only stains *Cryptosporidium* and *Giardia*. EasyStain™ does not stain typical contaminants of water (e.g. algae) which can be mistaken for C&G



EasyStain™ IgG1



Competitor IgM

Photographs courtesy of G.D. Di Giovanni and W.Q. Betancourt, TAES, Texas A&M University

Why use EasyStain™

- No background staining – reduces microscope time & increases laboratory throughput
- No false positives – does not stain algae, bacteria, yeast, proteins & minerals often present in water samples
- 15 minute incubation – rapid results
- Used by more laboratories worldwide than any other stain

Globally Approved

- United States Environmental Protection Agency (Method 1622/1623)
- United Kingdom Drinking Water Inspectorate
- Japanese Water Works Association

Water Industry Enquires:
BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com


A BIOMÉRIEUX COMPANY


BIOMÉRIEUX
INDUSTRY