

Batch Number: **B 112****Fluorescence analysis - Mean Fluorescence Intensity***Cryptosporidium* oocysts: 1226.61*Giardia* cysts: 1353.39**Non-specific binding assay**

Number of particles detected: 6

Control number of particles detected: 7

**General**

Safety: This product contains sodium azide. Handle with care.  
Storage: 2-8°C. Do not freeze.  
Catalogue Number: EST-CG20  
Preparation Date: 17 November 2015  
Expiration Date: 17 November 2016

**Quality Officer:**

Linh Le

**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

