

Batch Number: B 128**Fluorescence analysis - Mean Fluorescence Intensity***Cryptosporidium* oocysts: **1051.95***Giardia* cysts: **1217.61****Non-specific binding assay**Number of particles detected: **0**Control number of particles detected: **0****General**

Safety:	This product contains sodium azide. Handle with care.
Storage:	2-8°C. Do not freeze.
Catalogue Number:	EST-CG80-D, EST-CG20-D, W-EST-CG80-D, W-EST-CG20-D
Preparation Date:	5 September 2018
Expiration Date:	6 September 2019

Quality Assurance Specialist

Linh Le
13 September 2018

**Water Industry Enquires:**

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8221 9589
www.btfbio.com

