

Batch Number: **B 548**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 10 September 2015
Expiration Date: 8 January 2016
Volume: 0.591 ml +/- 27 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

	Mean	St.Dev.
Cryptosporidium count:	99	1.4

DAPI staining:

Cryptosporidium % +ve	100 %
------------------------------	--------------

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 22 August 2015
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Officer:



.....
Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

