

Batch Number: **B 581**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 10 November 2016
Expiration Date: 10 March 2017
Volume: 0.676 ml +/- 10 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

| Counts | | |
|-------------------------------|-----------|------------|
| <i>(Method Ref: ESCG-001)</i> | Mean | St.Dev. |
| Cryptosporidium count: | 99 | 1.7 |
| DAPI staining: | | |
| Cryptosporidium % +ve | 96 | % |

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 26 October 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

